REF ABV700

AbViser® IAP Monitor

OPERATIONS MANUAL

AbViser® IAP Monitor

Manufactured By:

Wolfe Tory Medical, Inc.

79 West 4500 South
Suite 18
Salt Lake City, Utah 84107
801-281-3000

www.wolfetory.com

EU Authorized Representative:

Cavendish Scott Ltd.
P.O. Box 107
SG5 1 FW, England
WWW.CSEAST.COM


© 2009 Wolfe Tory Medical, Inc. All rights reserved.

AbViser® is a registered trademark of Wolfe Tory Medical, Inc.
Made in the U.S.A.

30-0228 Rev. 11/09

CONTROL PANEL, CONNECTORS AND INDICATORS

(Refer to Monitor Illustrations on page 1):

#	DESCRIPTION
1	IAP Pressure Display [mmHg].
2	Display status: Interface cable connection, Record-Hold, and Zeroing.
3	Date [yyyy/mm/dd].
4	Time [24:00].
5	Remaining battery capacity.
6	Charging method: either supplied 5V DC medical grade power supply or USB connection.
7	RECORD - Freezes IAP reading on display and date/time stamps information.
8	- Turns the Monitor "ON" and "OFF".
9	MENU - Retrieves IAP History from memory and allows adjustment of internal real-time clock.
10	ZERO - "Zeros" the monitor when the pressure transducer is vented to atmospheric pressure. Clears the internally stored IAP history. Used to select adjustment field for adjusting internal clock.
11	- Scrolls through internally stored IAP History. Adjust internal real-time clock. Used as confirmation during clearing of internally stored IAP history.
12	Interface cable connector.
13	Power Connector.
14	USB Connector.
15	Mounting Bracket.
16	Pole-clamp.

SYMBOLS

SYMBOL	DEFINITION
	Attention, consult accompanying documents.
	Type BF Equipment.
	Signifies that this device has been tested and certified by Intertek to comply with applicable safety standards.
	CE 0050 indicates that the device meets the requirements of 93/42/EEC.
	Pressure interface cable and/or transducer not connected.
	Battery Charging.
	Battery Capacity indicator with multiple segments.
	Connected via USB.
	Connected to Power Supply.

Indications for Use

The AbViser® IAP Monitor is intended to display and record intra-abdominal pressure obtained from an AbViser® disposable monitoring device.

Contraindications

This monitor is not designed, sold, or intended for use except as indicated. Federal law (U.S.A.) restricts this device to sale by or on the order of a clinician.

Getting Started

Take a few minutes to ensure the following contents are included: monitor, power supply, power cable, USB cable, pole clamp, and operations manual.

Pressure transducer “interface-cables” are sold separately.

Before attempting to use the AbViser® IAP Monitor, become acquainted with the monitor and its accessories. Unpack the items carefully and check for any damage. If there is a problem contact Wolfe Tory Medical.

It is recommended that all the packaging be retained, in case the monitor must be returned for service in the future.

Power Sources

The AbViser® IAP Monitor is designed to operate from either the internal battery, or from the supplied power supply. Only use the AbViser® IAP Monitor with the supplied power supply.

Powering the Monitor

Before the AbViser® IAP Monitor is used for the first time, the battery should be charged for at least eight (8) hours. A “plug” icon is displayed when the monitor is connected to the power supply. A “USB” icon is displayed when the battery is connected to the PC via the USB connection. If the monitor is simultaneously connected to the power supply and the USB, only the “plug” icon shall be displayed. A “four-segmented-battery” icon indicates the remaining battery capacity. If the battery is flashing, the battery capacity is at a critically low state and the monitor shall turn off momentarily, it is strongly advised to recharge the battery to prevent any disruption.

Note: Connecting the power supply or USB connection during regular use does not interrupt patient IAP monitoring.

If the display does not show a “plug” icon when plugged into the power supply, check all connections and make sure the ac power source is on.

Using the Monitor

Attaching the Monitor from the Pole Clamp

1. The pole clamp must be secured to the pole, or any sturdy structure, with a thickness less than 1.5 [inches].
2. Align the bottom of the monitor mounting bracket to the lower arm of the pole clamp, as depicted in Figure 3.
3. With the lower arm of the pole clamp captured by the bottom of the mounting bracket, push the monitor in until a “click” is heard and the upper arm is captured by the upper mounting bracket.

Detaching the Monitor from the Pole clamp

1. Grasp the upper portion of the monitor and push down.
2. While still pushing down, pull the monitor outward, as depicted in Figure 5.

Turning the Monitor “ON” and “OFF”

To turn the monitor “ON”, touch the “” key. An audible confirmation with the pitch ascending low to high will sound when the monitor is turned “ON”. To turn the monitor “OFF”, touch the “” key. An audible confirmation will descend from high to low when the monitor is turned “OFF”.

Internal Real-Time Clock Adjustment

The AbViser® IAP Monitor contains an internal real-time clock. If the AbViser® IAP Monitor is connected to a computer via the USB interface, the internal real-time clock can synchronize with the connected computer’s real-time clock. The clock may also be manually adjusted via the key pad.

Internal real-time clock adjustment procedure.

1. Touch the “MENU” key twice to enter the manual clock adjustment mode. An audible confirmation will sound with one (1) click each time the “MENU” key is touched.
2. Touch the “ZERO” key to toggle to the next adjustment field. The first field selected for adjustment is the hour. An audible confirmation will sound with one (1) beep when the “ZERO” key is touched.
3. Touch the “” key to adjust the highlighted adjustment field. An audible confirmation will sound with one (1) click when the “” key is touched.
 - a. Individual touches will incrementally adjust the selected field up.
4. Touch the “RECORD” key to save the current clock settings or touch the “MENU” key to discard the current clock settings and return to the home screen.

Obtaining Intra-Abdominal Pressure (IAP) on Display

The monitor must be turned “ON” and zeroed before an accurate IAP measurement can be made.

1. Inspect cables and connectors for damage; replace the cables when aging, tearing, or a loose connection is apparent. Replacement cables can be obtained from Wolfe Tory Medical.
2. Connect the “interface-cable” to the IAP connector on the monitor.
3. Connect the transducer to the “interface-cable.”
4. Follow instructions per AbViser® IAP disposable device for zeroing the transducer. Once ready for zeroing, touch the “ZERO” key on the monitor.
5. Follow instructions per AbViser® IAP disposable device for measuring IAP.
6. The intra-abdominal pressure will appear on the display in [mmHg].
7. Touch the “RECORD” key to record the current IAP measurement and date/time stamp information.

Note: If the “RECORD” key is not touched after three (3) minutes elapses from the moment the infusion bolus is introduced, or the “ZERO” is established, the monitor will automatically display: “READY FOR NEW IAP MEASUREMENT” and does not need to be rezeroed.

Record Data

The AbViser® IAP Monitor has the ability to record the current IAP reading on the display and store the information internally for twenty-four (24) measurements.

1. The recorded information shall be the IAP pressure [mmHg], date [yyyy/mm/dd] and time [24:00] at the moment the “RECORD” key is touched.

2. This reading shall remain displayed on the monitor in the "IAP RECORDED-HOLDING" mode until the monitor detects an infusion pressure spike from subsequent IAP testing, or when the operator touches the "RECORD" key.
3. Once twenty-four measurements are recorded in the internal memory and a new measurement is recorded, the oldest measurement shall be overwritten by the latest measurement.

Reviewing Internally Stored Data

The AbViser® IAP Monitor allows the user to review the twenty-four (24) internally stored IAP measurements.

1. Touch the "MENU" key once to recall the last twenty-four recorded IAP measurements and time stamp information. The data is presented in tabular form, with each measurement residing on a given row. The first column is the date stamp information, the second column is the time, and the third column is the IAP measurement.
2. To scroll through the table one event at a time, touch the "→" key. The table shall return to the beginning after the 24th data measurement has been displayed.
3. If the "MENU" key is touched, the monitor will return to the home screen. Also, if no activity is detected within thirty (30) seconds, the monitor will automatically return to the home screen.

Note: The monitor will automatically display the following message after three (3) minutes elapses from the moment the monitor exits the IAP History and the monitor has been previously "Zeroed":
"READY FOR NEW IAP MEASUREMENT."

Erasing Internally Stored Data

1. Select the "MENU" key to review the internally stored data.
2. Touch the "ZERO" key. A confirmation screen "CLEAR DATA LOG? →=Y" shall be displayed.
3. Touch the "→" key as confirmation to erase all of the IAP history from the internal memory, or touch "MENU" key to exit without erasing the stored IAP History.

CAUTION: Always erase the IAP history between patients.

Establishing connection to a Computer

The appropriate driver must be installed before the AbViser® IAP Monitor is able to transfer all the information stored to a computer via the USB interface.

1. Download the ABV700 driver from the following URL:
http://wolfetory.com/abviser_monitor/downloads/
2. Turn the IAP monitor "ON."
3. Connect the USB data cable to the USB socket on the monitor; ensure that the USB terminator interfaces correctly with the USB socket. Refer to the illustration of the bottom of the monitor.
4. Connect the other end of the USB data cable to the computer.
5. Follow Hospital protocol/policy for installing drivers.

Transferring Internally Stored Data to a Computer

The AbViser® IAP Monitor has the ability to transfer all the information stored internally to a computer via the USB.

1. Download and install ABV700 data acquisition software, available at the following URL:
http://wolfetory.com/abviser_monitor/downloads/
2. Run the ABV700 software.

Maintenance

1. Cleaning the monitor and cable inspection: Before cleaning, thoroughly inspect the monitor and all accessories for any signs of damage, or cracks. While gently bending and flexing cables, inspect for damage such as cracks, cuts, abrasions, extreme wear, exposed wires or bent connectors. Confirm connectors securely engage. Replacement cables may be acquired from Wolfe Tory Medical. The monitor is not end user serviceable, contact Wolfe Tory Medical for repair service.
 - a. The exterior of the monitor may be wiped clean with a soft cloth slightly dampened with water. Thoroughly wipe off any excess water. Do not let water run into connector openings or crevices. The monitor and accessories should be periodically inspected for wear and damage.
 - b. Do NOT immerse the unit.
 - c. Do NOT clean with Isopropyl alcohol, or other solvents.
 - d. Do NOT immerse the accessories.
2. Storage
 - a. If it becomes necessary to store the monitor for an extended period of time, first fully charge the battery. Store the monitor and the power supply in the original packaging materials.
 - i. Batteries should always be fully charged before being placed in storage. It is best to keep the battery charged as fully as practical and never store the monitor with the battery in a discharged condition. Battery charging will take place as long as the monitor remains connected to an external AC power source.
 - ii. To ensure that the battery will be ready for portable operation, keep the monitor connected to the power supply whenever possible.
 - iii. Repeated failure to fully charge the battery will result in a significant reduction in battery life.
3. Batteries
 - a. Batteries are NOT user replaceable. Contact Wolfe Tory Medical for service.
 - b. Repeated failure to fully charge the battery will result in a significant reduction in battery life.
4. Disposal
 - a. When the device reaches end-of-life, follow hospital disposal protocol for electronic devices containing rechargeable batteries.

Warranty & Servicing

Wolfe Tory Medical, Inc. ("WTM") warrants to the original end user customer that its products are free from defects in material and workmanship. Subject to the conditions and limitations set forth below, WTM will, at its option, either repair or replace any part of its products that prove defective by reason of improper workmanship or materials. Repaired parts or replacement products will be provided by WTM and will be either new or refurbished to be functionally equivalent to new.

This limited warranty does not cover any damage to this product that results from improper installation, accident, abuse, misuse, natural disaster, insufficient or excessive electrical supply, abnormal mechanical or environmental conditions, or any unauthorized disassembly, repair, or modification. This limited warranty also does not apply to any product on which the original identification information has been altered, obliterated or removed, has not been handled or packaged correctly, has been sold as second-hand or has been resold contrary to the US export regulations.

This limited warranty covers only repair, or replacement for defective WTM products, as

provided above. WTM is not liable for, and does not cover under warranty, any loss of data or any costs associated with determining the source of system problems or removing, servicing or installing WTM products. This warranty excludes 3rd party software, connected equipment or stored data. In the event of a claim, WTM's sole obligation shall be replacement of the hardware.

Duration of Warranty

The monitor is covered under a one-year warranty from the date of purchase. This limited warranty applies only to the original end user customer of the product for so long as the original end user customer owns the product. This limited warranty is non-transferable.

Return Material Authorization (RMA) Policy

1. Customer will assume all shipping charges on items returned to WTM. It is recommended to ship merchandise with insurance or to send by any carrier that provides tracking/delivery confirmation/insurance.
2. WTM is not responsible for any lost, or damage of the merchandise during/caused by the shipment.
3. WTM will assume all shipping charges on items returned to customers if service is covered by the warranty; otherwise customer is responsible for all shipping charges for items returned to the customer.
4. Defective items will be repaired or replaced at no charge provided that items are not abused, misused, neglected or modified.
5. WTM will promptly repair or replace items within thirty (30) business days.
6. A valid RMA number must accompany any merchandise returned to WTM. An RMA number may be obtained by contacting WTM customer service at 1-801-281-3000.

Product Compliance

The AbViser® IAP Monitor is classified in the following categories for compliance with EN 60601-1:

- Internally powered by batteries
- Transportable
- For continuous operation
- Class 1 Medical Grade Power Supply
- Type BF Device (Body Floating type device that will be electrically connected to the patient via cables or liquid filled catheters.)

Monitor IAP Limits of Accuracy Specification

Characteristic	Specification
IAP Range	-30 mmHg to 99 mmHg
IAP Accuracy	±1 mmHg plus Transducer Error (0 -50 mmHg) ±5 mmHg plus Transducer Error (51- 99 mmHg)

Monitor Specification

Characteristic	Specification
Type of protection against electric shock-Power Supply	Power supply class 1
Type of protection against electric shock-Monitor	Protective earth not available in monitor.
Method of disinfection	Not suitable for autoclaving, see cleaning instructions
Presence of Flammable Mixtures	Not suitable for use in presence of flammable mixtures

Battery Specification

Characteristic	Specification
Mode of operation	Continuous
Battery Pack Type	Lithium Ion
Battery Pack Capacity	4400 mAh
Battery Recharger circuitry	Internal, powered by external power supply
DC Input power required	3.7V DC
Typical Battery Life	48 hours

Monitor Environmental Specifications

Characteristic	Specification
Operating Temperature	0° to 40° C
Shipping & Storage Temperature	-20° to 60° C
Operating Altitude	0 to 15,000 ft (0 to 4,572 m)
Shipping and Storage Altitude	0 to 40,000 ft (0 to 12,192 m)
Operating Relative Humidity	15% to 95%, noncondensing per MIL STD 810E, Procedure 1-natural
Shipping and Storage Relative Humidity	15% to 95%, noncondensing per MIL STD 810E, Procedure 1-natural

Power Supply Specification

Characteristic	Specification
Type of protection against electric shock	Class 1 Medical Grade External Power Supply
Degree of protection against harmful ingress of water	For ordinary, indoor locations only.
Method of disinfection	Not suitable for autoclaving
Presence of Flammable Mixtures	Not suitable for use in presence of flammable mixtures
Operating Temperature	0° to 40° C with no de-rating
Shipping & Storage Temperature	-30° to 85° C
Operating Altitude	0 to 10,000 ft (0 to 3,048 m)
Shipping & Storage Altitude	0 to 40,000 ft (0 to 12,192 m)
Operating Relative Humidity	5% to 95%, noncondensing
Shipping & Storage Relative Humidity	5% to 95%, noncondensing
Rated Input	100-240 VAC ± 10%, 47-63 Hz
Rated Output	5VDC, 2.0A
Additional Features	Detachable power cord

Warnings & Disclaimers:

- The IAP monitor is to be operated only by qualified personnel. The operator of this IAP monitor should read this entire manual prior to initial use.
- If the accuracy of any reading is questionable, first check for the proper setup of the coupled AbViser® IAP device per the AbViser® IAP Monitoring Kit Instructions. If needed, check the AbViser® IAP Monitor for proper setup per operating manual.
- Motion artifact can affect the accuracy of the IAP measurement. Minimize patient motion whenever possible.
- Do not use the monitor in an MRI suite or a hyperbaric chamber, or expose to X-ray energy.
- Do not sterilize the monitor or accessories, damage will occur by sterilization.
- Only use this monitor with Wolfe Tory Medical AbViser® Intra-abdominal Pressure (IAP) Disposable monitoring devices.
- Place the monitor and accessories in locations where they cannot harm the patient if they should fall from their shelf or mount.
- Carefully route the patient cabling to reduce the possibility of patient entanglement or strangulation.
- To ensure patient safety, accessories must be used according to your hospital's protocols and the manufacturer's recommendations. Always refer to the manufacturer's directions for use.
- Only use the monitor with the supplied power supply. Inspect the power supply cord periodically for fraying or other damage, and replace the power supply as needed. Contact WTM for replacement power supply. Do not operate the monitor with a damaged power supply cord or plug.
- When using the power supply with this monitor, be sure to connect the power supply to a three-wire grounded, hospital-grade receptacle (North America only). Do not under any circumstances attempt to remove the grounding conductor from the power plug of the power supply. Do not plug the power supply into an extension cord. If there is any doubt about the integrity of the protective earth ground of the receptacle for the power supply, contact your engineering department for assistance in identifying the proper power receptacle and making appropriate power connections.
- Make frequent electrical and visual checks on all cables.
- Do not lift the monitor by any attached cables.
- The monitor does not include any user-serviceable components. Do not attempt to open the monitor to perform any service on the monitor.
- A product that has been dropped or severely abused should be checked by qualified service personnel to verify proper operations. Contact WTM for repair or replacement.
- A Foley Catheter is not supplied with the monitor. Foley catheters are supplied by the hospital. Foley catheters are single-patient use devices (Ⓢ) and are EtO sterilized STERILE EO. Please follow the instructions provided by the Foley Catheter manufacturer and medical facility guidelines to ensure patient safety and product sterility.
- Wolfe Tory Medical assumes no responsibility for usage not in accordance with this manual that results in illegal or improper use of the monitor.
- Wolfe Tory Medical will not be responsible for any injury to the user or other person(s) that may result from accidents during use of the monitor.
- No implied license: possession or purchase of this device does not convey any expressed or implied license to use the device with unauthorized accessories.
- This manual may be wholly or partially subject to change without notice.
- All rights are reserved. No one is permitted to reproduce or duplicate the whole or any part of this manual without written consent from Wolfe Tory Medical.

- In no case will WTM be liable to the purchaser, or to any user of the WTM product, for any losses including, expenses, revenues, profits, equipment, or any other subsequent damage that might arise from the use or misuse of any WTM product, which also applies in the situation where WTM has been previously warned of any likely adverse effects occurring that might cause any of the above situations to transpire.
- This disclaimer as stated above supersedes any other expressed or implied warranties whether it is written or oral. Except as expressed above, WTM expressly disclaims all warranties not stated herein, including, to the extent permitted by applicable law, any implied warranty of merchantability or fitness for a particular purpose.

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